Startup Sequence

System User Manual

Addendum

Intuitive Surgical, Inc.
1266 Kifer Road
Sunnyvale, CA 94086 USA
intuitive.com

Intuitive Surgical SAS
11 avenue de Canteranne
33600 Pessac, France

PN 557231-01 Rev. B 2024.03
(Multilingual)

Copyright

© 2024 Intuitive Surgical Operations, Inc. All rights reserved.

Trademarks

Product and brand names/logos are trademarks or registered trademarks of Intuitive Surgical or their respective owner. See intuitive.com/trademarks.

Rx only

Contents

English (en)	4
Česky (cs)	5
Dansk (da)	6
Deutsch (de)	7
Español (es)	8
Français (fr)	10
Hrvatski (hr)	12
Italiano (it)	13
Magyar (hu)	14
Nederlands (nl)	15
Norsk (no)	16
Polski (pl)	17
Português (pt)	19
Português Brasileiro (pt-BR)	21
Română (ro)	22
Slovenčina (sk)	24
Slovenščina (sl)	25
Suomi (fi)	26
Svenska (sv)	27
Tiếng Việt (vi)	28
Türkçe (tr)	29
Ελληνικά (el)	30
Български (bg)	32
Русский (ru)	34
한국어 (ko)	35
日本語 (ja)	36
繁體中文 (zh-TW)	37

English (en)

System User Manual Addendum

This addendum to the Da Vinci Xi System User Manual and the Da Vinci X System User Manual contains an update to the system startup sequence information found in Chapter 5. For full system information, including General Information and Professional Instructions for Use, see the applicable system user manual. Keep this addendum with your system documentation.

5.2 Startup Sequence

During the startup sequence, the system performs an integrity test. The Patient Cart arms and Surgeon Console hand controls perform various movements during this integrity test. An audible beep and voice annunciated feedback (if enabled) indicates that the system is ready.

Patient Cart

During system startup, the Patient Cart arms move and perform a short mechanical integrity test. During the test, the arm carriage LEDs illuminate pulsing blue. Once the system integrity test is successfully completed, all arm LEDs illuminate solid blue, and the Power button LEDs illuminate solid blue when powered on.

- If the system detects an instrument or cannula before startup, the arms will not move during startup. This feature allows the system to stay connected to a patient during the startup sequence.
- If an arm bumps into something during the test, use the instrument or port clutch buttons to move the arm clear of the obstacle.

Surgeon Console

The hand controls perform a self-test. The hand controls move to their start position and must arrive there for the system to work.

- If a hand control is impeded in some way, move it by hand to free it and it will move to its start position.
- During the startup sequence, do not put anything in the 3D viewer.
- Do not activate any of the system controls, including the clutch buttons, foot pedals, etc.
 - While most button activations will be ignored during the startup sequence, some may cause a non-recoverable fault, which means that the system will have to be restarted. See Non-recoverable Faults in Appendix B for more information.
- The Emergency Stop feature is available during startup if needed. Once the system beeps, the surgeon may interact with the system controls. See Emergency Stop – Surgeon Console in Appendix B for more information.

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