




## What's new in Arthroscopy and Sports Medicine?

### DrillCut-X® ARTHRO



The new DRILLCUT-X® ARTHRO shaver handpiece stands out due to its excellent cleanability and exceptional versatility. The distal end of the shaver handpiece is detachable so that optimal cleaning of previously hard-to-access areas is possible. Another striking feature of the DRILLCUT-X® ARTHRO shaver handpiece is its ergonomic design. Due to its small dimensions, the tool is suitable for use in both large and small joints. The shaver handpiece is controlled via the UNIDRIVE® S III ARTHRO console.

**Product Features:**

- Optimal reprocessing of previously hard-to-reach areas thanks to disassembly feature
- Sophisticated ergonomic design allows various hand positions
- Possible to lock blades in 4 positions (0°, 90°, 180°, 270°)
- Suction regulation via a rocker lever
- Low-noise operation
- Handpiece weighs 310 g



Optimized reprocessing thanks to unique disassembly feature

**Order information:**

28200 DX      **DRILLCUT-X® ARTHRO Shaver Handpiece**, up to 8000 rpm, for use with UNIDRIVE® S III ARTHRO SCB, as of Software Status 1.10.

**It is recommended to check the suitability of the product for the intended procedure prior to use.**

*For further information about the mentioned products, please contact  
KARL STORZ GmbH & Co. KG, Mittelstraße 8, 78532 Tuttlingen/Germany  
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105, E-Mail: [info@karlstorz.com](mailto:info@karlstorz.com)  
[www.karlstorz.com](http://www.karlstorz.com)*