



What's new in Anesthesiology and Emergency Medicine?

LED BRITE LITE Battery Light Source

Accessories for intubation fiberscopes and endoscopes



Special Features:

- Battery light source with extremely high light intensity > 100 lm / > 150 klx
- Available as battery and rechargeable version
- Absolute white light due to LED technology
- Special light focus allows optimal light adjustment at the endoscope connector
- LED service life up to 50,000 hours
- Burning time over 120 min
- Waterproof, fully immersible for cleaning and disinfection (11301 D1/D3)



11301 D1/D3/DE/DF



11301 DG

- 11301 D1 **Battery Light Source LED for Endoscopes**, with fine screw thread, brightness > 100 lm / > 150 klx, burning time > 120 min, weight approx. 150 g, waterproof and fully immersible for manual cleaning and disinfection, with 2 Photo Batteries 121306 P
- 11301 D3 **Same**, with coarse thread
- 121306 P **Photo Battery**, lithium, 3 V, CR 123 A
- 11301 DE **Battery Light Source LED for Endoscopes**, rechargeable, with click connection, brightness > 110 lm / > 150 klx, color temperature 5500 K, lithium-ion batteries, charging time 60 min, burning time at 100% brightness 40 min, weight approx. 150 g, suitable for wipe disinfection
- 11301 DF **Same**, with fast screw thread
- 11301 DG **Charging Unit**, for 11301 DE/11301 DF, for two LED battery light sources, with fixed integrated power supply and adaptor for EU, UK, USA and Australia, power supply 110-240 VAC, 50/60 Hz, suitable for surface disinfection
- 11301 DH **Holder**, for Charging Units 11301 DG, 8546 LE and 8401 XDL

It is recommended to check the suitability of the product for the intended procedure prior to use.

Consent to receive electronic information

Yes, I agree to receive future information by email at the following address:

_____	_____
Email	Name
_____	_____
Department / Practice	Street address
_____	_____
ZIP, Town	Signature

I agree to my data being stored at KARL STORZ for this purpose. I can withdraw my consent at any time and without giving reasons by emailing KARL STORZ at info@karlstorz.com. KARL STORZ will not make these data available to third parties.

*For further information about the mentioned products, please contact
KARL STORZ GmbH & Co. KG, Mittelstraße 8, 78532 Tuttlingen/Germany
Phone: +49 (0)7461 708-0, Fax: +49 (0)7461 708-105, E-Mail: info@karlstorz.com
www.karlstorz.com*